

Q4₂



Q4₂
shown in *Red Corvette* with TB3 seating
system and Q-Logic 3e Advanced Drive
Controls with Bluetooth®

With features like standard 6.25 mph motors and advanced responsiveness, the Q4₂ is engineered to meet the performance needs of the most active user.

The Q4₂ accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value.

QUANTUM®

Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic 14" Knobby
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
MAXIMUM SPEED ¹	Up to 6.25 mph / 10 km/h
GROUND CLEARANCE	3.03" (at battery tray)
TURNING RADIUS ^{2,4}	20.5"
OVERALL LENGTH ^{2,4}	35.5" (33.74" Stowage Length)
BASE WIDTH ⁴	24.41"
BATTERY SIZE	NF-22
RANGE ^{1,3}	Up to 19.2 mi / 30.8 km
BATTERY WEIGHT ⁵	38 lb / 17.24 kg
STANDARD ELECTRONICS	QLNES
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Motor Mid-Wheel 6 [®]
SUSPENSION	ATX (Active-Trac [®] with extra stability)
WEIGHT CAPACITY	300 lb / 136 kg
BASE WEIGHT	138.6 lb / 63 kg
SEATING SIZES ⁴	
TRU-Balance[®] 3 Static Seat	Width: 12-22" Depth: 12-22"
TRU-Balance[®] 3 Power Tilt	Width: 12-22" Depth: 12-22"
Q-Captain's Seating	Width: 16-20" Depth: 16-22"
TRU-Balance[®] Flex Seat	Width: 15-22" Depth: 15-20"
SEAT-TO-FLOOR HEIGHTS ⁴	
TRU-Balance[®] 3 Static Seat	16.5-18.5"
TRU-Balance[®] 3 Power Tilt ³	17.5-19.5"
TRU-Balance[®] Flex Seat	16.5-18.25"

1. Varies with user weight, terrain type, battery charge, battery condition and tire pressure.
2. Without front riggings
3. Battery capacity may change when adding power positioning systems.
4. All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
5. Battery weight may vary based on manufacturer.


Features

- Standard 6.25 mph motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Q-Logic QLNES drive control with Bluetooth connectivity
- ATX Suspension
- Mid-wheel drive stability
- Bold aesthetics

Options

- TRU-Balance[®] 3 seating
- Q-Captain's seating
- Flex seating
- Stealth Seating and Positioning
- Vent trays
- 4" wide knobby extended traction tires
- LED fender lights
- Coloured wheel inserts
- Upgradeable Q-Logic and Stealth drive controls
- Accu-Trac Advanced Tracking Technology
- Available with WC19 crash-test approved occupied transit kit option on tilt & recline and iLevel
- WC19 approved transit kit.

Colors



QUANTUM[®]

#1 FOR REHAB POWER

quantumrehab.ca

