

# QUANTUM® R-TRAK

REAR-WHEEL DRIVE POWER CHAIR

**R-TRAK**  
shown in *Evolution Orange* with  
TRU-Balance® 3 power positioning,  
Q-Logic 3 drive controls and  
Stealth® TRU-Comfort 2 cushion.



## NEXT LEVEL REAR-WHEEL DRIVE

The R-TRAK revolutionizes the rear-wheel drive power chair experience. Independent SRS suspension on all four wheels combined with a compact base provide a smoother, more maneuverable driving experience. Standard features include fender lights, USB charger and drive wheel color accents. The R-TRAK is also available with our 10" safe seat elevation, allowing the client to drive up to 3.5 mph while elevated.



**QUANTUM®**  
#1 FOR REHAB POWER



| [quantumrehab.ca](http://quantumrehab.ca) |

## Specifications

<b>DRIVE WHEELS</b>	Standard: 14" pneumatic Optional: 14" flat free 14" knobby
<b>CASTER WHEELS</b>	Front: 9" flat free
<b>ANTI-TIP WHEELS</b>	Rear: 4" solid
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 11 km/h (7 mph) (3.5 mph when elevated)
<b>GROUND CLEARANCE</b>	3.97" (battery tray)
<b>TURNING RADIUS</b> <sup>2</sup>	25.7"
<b>OVERALL LENGTH</b> <sup>2</sup>	34.8"
<b>BASE WIDTH</b>	24.2"
<b>SEATING SIZES</b>	
<b>TRU-Balance® 3 Static Seat</b>	Width: 12-22" Depth: 12-22"
<b>TRU-Balance® 3 Power Tilt</b>	Width: 12-22" Depth: 12-22"
<b>TRU-Balance® 3 Power Lift &amp; Tilt</b>	Width: 12-22" Depth: 12-22"
<b>TRU-Balance® Flex Seating</b>	Width: 15-22" Depth: 15-20"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>TRU-Balance® 3 Static Seat</b> <sup>5</sup>	16-18"
<b>TRU-Balance® 3 Power Tilt</b> <sup>5</sup>	18.5-19.5"
<b>TRU-Balance® 3 Power Lift &amp; Tilt</b> <sup>5</sup>	17.25-18.25"
<b>TRU-Balance® Flex Seating</b>	16.5-17.5"
<b>BATTERY SIZE</b> <sup>3</sup>	Group 24 (2)
<b>BATTERY WEIGHT</b> <sup>7</sup>	Group 24: 23.4 kg (51.5 lb) each
<b>AVAILABLE ELECTRONICS</b>	120A Q-Logic 3 NE+
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	4-pole motors
<b>WEIGHT CAPACITY</b> <sup>4</sup>	136 kg (300 lb)
<b>BASE WEIGHT</b>	58 kg (128 lb)
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard: (K0848 <sup>^^</sup> /K0849 <sup>^^</sup> ) Single Power: (K0856 <sup>^^</sup> /K0857 <sup>^^</sup> ) Multiple Power: (K0861 <sup>^^</sup> )

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 7) All Quantum Rehab® units can be ordered with specialty controls.
- 8) Call for availability
- 9) Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front 2's suspension system, seat-to-floor heights are approximately 0.25" lower when the unit is loaded.

<sup>^^</sup>Pending PDAC code verification.

## Standard Features

- 4-pole motors
- Maximum speed of up to 11 km/h (7 mph)
- LED lighting in front caster arms and fenders
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions


## Options

- 10" of power adjustable safe seat elevation at speeds of up to 5.6 km/h (3.5 mph)
- Quantum® backup camera
- Trilogy vent tray <sup>8</sup>
- Compact vent tray <sup>8</sup>

## Colors



## Hydrographics



**QUANTUM®**

#1 FOR REHAB POWER

quantumrehab.ca

