



AVAILABLE WITH 12" iLEVEL® AT 3.5 MPH



AVAILABLE WITH 12" iLEVEL® AT 4.5 MPH



AVAILABLE WITH 10" ELEVATION® AT 3.5 MPH



AVAILABLE WITH 10" ELEVATION AT 3.5 MPH



AVAILABLE WITH 10" HD iLEVEL® AT 3.5 MPH



AVAILABLE WITH 10" ELEVATION AT 2.5 MPH



	J6	Q4	EDGE 3 STRETTO®	EDGE 3	R-TRAK	4FRONT® 2	EDGE 9 HD	4FRONT® 2 HD	Q1450
	Powerr BaseType 1: Adult: WP0000564	Powerr BaseType 2: Adult: WP0001541	Power Base Type 3: Adult: WP0001806 Paediatric: WP0001807	Powerr BaseType 3: Adult: WP0001575 Paediatric: WP0001576	Powerr BaseType 3: Adult: WP0001987 Paediatric: WP0001988	Powerr BaseType 3: Adult: WP0001441 Paediatric: WP0001442	Powerr BaseType 3: Adult: WP0001058 Paediatric: WP0001059	Powerr BaseType 3: Adult: WP0001900	Powerr BaseType 3: Adult: WP0001021
DRIVE WHEEL CONFIGURATION	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Rear-wheel Drive	Front-wheel Drive	Mid-wheel Drive	Front-wheel Drive	Front-wheel Drive
WEIGHT CAPACITY ⁷	300 lb / 136 kg ⁹	300 lb / 136 kg	300 lb / 136 kg	300 lb / 136 kg	300 lb / 136 kg	300 lb / 136 kg	450 lb / 204 kg ³	450 lb / 204 kg	600 lb / 272 kg
MAXIMUM SPEED ⁴	Up to 4 mph / 6.4 km/h	Up to 5 mph / 8 km/h Optional: Up to 6 mph / 9.6 km/h ^{11,13}	Up to 6.25 mph / 10 km/h ¹³ (3.5 mph / 5.6 km/h with iLevel®)	Up to 6.25 mph / 10 km/h ¹³ (4.5 mph / 7.2 km/h with iLevel®) Optional: Up to 7 mph / 11.27 km/h ¹	Up to 7 mph / 11.27 km/h (3.5 mph / 5.6 km/h when elevated)	Up to 6.25 mph / 10 km/h (3.5 mph / 5.6 km/h when elevated)	Up to 4.7 mph / 7.5 km/h (3.5 mph / 5.6 km/h with iLevel®) ³	Up to 5 mph / 8 km/h (2.7 mph / 4.3 km/h when elevated)	Up to 4.5 mph / 7.2 km/h
GROUND CLEARANCE ⁸	2" (caster arm)	2.9" (battery tray)	2.9"	2.7" (battery tray)	3.97" (battery tray)	3.14" (battery tray)	3" (base frame)	3.14" (battery tray)	1.5" (anti-tip wheels)
TURNING RADIUS ^{5,8}	22"	20.5"	21"	20.5"	25.7"	25.9"	19.8"	25.9"	26.5"
LENGTH ^{5,8}	34.7"	35.5"	32.5"	35.5"	34.8"	39.9"	35.6"	39.9"	38.3"
WIDTH ⁸	23.2"	24.5"	21.5" (14" drive wheels) 20.5" (12.5" drive wheels)	24.4"	24.2"	24"	26.3"	24"	29.3"
DRIVE WHEELS	10" Flat Free Optional: 10" Pneumatic	14" Flat Free Optional: 14" Pneumatic	14" or 12.5" Pneumatic or Flat Free tires Optional: 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Flat Free
CASTER WHEELS	Front: 5" Solid Rear: 6" Solid	Front: 6" Solid Rear: 6" Solid	Front: 5" Solid Rear: 6" Solid	Front: 6" Solid Rear: 6" Solid	Front: 9" Flat Free	Rear: 9" Flat Free	Front: 6" Solid Rear: 6" Solid	Rear: 9" Flat Free	Rear: 9" Flat Free
ANTI-TIP WHEELS	N/A	N/A	N/A	N/A	Rear: 4" Solid	Front: 4" Solid	N/A	Front: 4" Solid	Front: 4" Solid
BASE WEIGHT	100.2 lb / 45.4 kg	131.6 lb / 59.7 kg	126.4 lb / 57.3 kg	138.6 lb / 62.8 kg	128 lb / 58.1 kg	130 lb / 59 kg	143.5 lb / 65 kg	130 lb / 59 kg	195 lb / 88.5 kg
BATTERY WEIGHT ¹⁴	23.4 lb / 10.6 kg each	40 lb / 23.6 kg each	40 AH: 31 lb / 14.1 kg each NF-22: 40 lb / 23.6 kg each	40 lb / 23.6 kg each	51.5 lb / 23.4 kg each	Group 24: 51.5 lb / 23.4 kg each Group 34: 44.8 lb / 20.3 kg each	51.5 lb / 23.4 kg each	Group 24: 51.5 lb / 23.4 kg each	51.5 lb / 23.4 kg each
MOTOR PACKAGES	2-pole	2-pole	2-pole	Standard: 2-pole Optional: 4-pole	4-pole	4-pole	4-pole	4-pole	4-pole
BATTERIES ⁶	U-1 (2)	NF-22 (2)	40AH (2) or NF-22 (2)	NF-22 (2)	Group 24 (2)	Group 24 (2) or Group 34 (2)	Group 24 (2)	Group 24 (2)	Group 24 (2)
BATTERY CHARGER	5A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board
AVAILABLE ELECTRONICS ¹²	Q-Logic NE Q-Logic NE+	Q-Logic 3EX Q-Logic 3e Q-Logic NE+	Q-Logic 3EX Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e QLNES	Q-Logic 3EX Q-Logic 3e Q-Logic NE Q-Logic NE+	Q-Logic 3EX Q-Logic 3e	Q-Logic 3EX Q-Logic 3e VR2 6-key
SEATING SIZES									
TRU-BALANCE® 4	N/A	N/A	W: 16-22" D: 14-22"	W: 16-22" D: 14-22"	N/A	W: 16-22" D: 14-22"	N/A	N/A	N/A
SYNERGY® HD STATIC	N/A	N/A	N/A	N/A	N/A	N/A	W: 20-28" D: 15-28"	N/A	W: 20-32" D: 20-28"
TRU-BALANCE® 3 STATIC	W: 14-20" D: 14-20"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 20-32" D: 20-25"	W: 20-24" D: 20-24"	N/A
TRU-BALANCE® 3 POWER TILT	W: 14-20" D: 14-20" ⁹	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 20-32" D: 20-25"	W: 20-24" D: 20-24"	N/A
TRU-BALANCE® 3 LIFT & TILT	N/A	N/A	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 20-24" D: 20-25" ¹⁰	W: 20-24" D: 20-24"	N/A
TRU-BALANCE® BARIATRIC POWER TILT	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	W: 20-32" D: 20-28"
TRU-BALANCE® FLEX SEAT	N/A	W: 15-22" D: 15-20"	W: 15-22" D: 15-20"	W: 15-22" D: 15-20"	W: 15-22" D: 15-20"	W: 15-22" D: 15-20"	N/A	N/A	N/A
SEAT-TO-FLOOR HEIGHTS ²									
TRU-BALANCE® 4 ⁸	N/A	N/A	18-19"	17.5-19.25"	N/A	17.75-18.75"	N/A	N/A	N/A
STATIC ⁸	16.75"	16.5 - 18.5"	16-18"	16.5 - 18.5"	16-18"	15 - 17"	16.5 - 18.5"	15-17"	16.5-18.5"
POWER TILT ⁸	16.75" ⁹	N/A	17.5-19.5"	17.5 - 19.5"	18.5-19.5"	17.75-19.75"	18.5-20.5"	17.75-19.75"	19.5 - 20.5"
LIFT & TILT ⁸	N/A	N/A	18-20"	17.5 - 19.5"	17.25-18.25"	17.75-19.75"	18.5 - 20.5"	17.75-19.75"	N/A
TRU-BALANCE® FLEX SEAT ⁸	N/A	16.5 - 18.25	16.5-18.5"	16.5" - 18.25"	16.5-17.5"	16.5 - 18.5"	N/A	N/A	N/A

FOOTNOTES:
 1. 7 mph can only be achieved with optional 4-pole motors and Q-Logic 3.
 2. Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front's suspension system, seat-to-floor heights are approximately 0.25" lower when the unit is loaded.
 3. Reduced to 400 lbs. with iLevel® power adjustable seat height.
 4. Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
 5. Without front riggings

6. Battery capacity may change when adding power positioning systems.
 7. Weight capacity is dictated by the selection of seating.
 8. All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Depths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
 9. 250 lbs. weight capacity with TRU-Balance® Systems

10. A minimum seating size of 20" x 20" is required with iLevel® on the Q6 Edge® HD; maximum seating size is 24" x 24"
 11. Standard with Accu-Trac and/or Q-Logic 3
 12. All Quantum units can be ordered with specialty controls
 13. Available only with Q-Logic 3.
 14. Battery weight may vary +/- 3 lbs. based on manufacturer.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

