



	TRU-Comfort 2 SPP	Essence SPP / Essence Plus SPP	Solution® SPP	Spectrum® Foam/Gel SPP	Premiere P	Glacial SP	Zen SP	Simplicity G
DESCRIPTION¹	 Multi-layered foam Deep contour base, CoolCore® cushion cover (has "pockets" on sides)	 Low profile contoured molded foam base, lateral and pre-ischial shelf, coccyx and sacral relief, CoolCore® cushion cover (has "pouch" in front), standard ABS rigidizer Essence Plus has adjustable 1/2" & 1" positioning wedges standard and optional FRG topper	 Moderate contour base foam, twin cell aqueous viscous gel insert in ischial well, CoolCore® cushion cover	 Low contour base foam with sacral relief, blue XXSoft FRG foam topper (sizes 22 wide and wider, green MedSoft FRG foam topper), CoolCore® cushion cover, urethane incontinence cover Spectrum® Foam with multi layer Visco foam insert in ischial well Spectrum® Gel with twin cell ultra-lightweight viscous gel insert in ischial well	 Moderate contour base foam with high density positioning wedges built in for offloading, gray XSoft FRG foam topper, CoolCore® cushion cover	 New SIM G base, Stealth Tek, 2-way stretch breathable cover (fluid resistant top and side panels), TPU heat-sealed low profile lightweight cooling gel topper	 Low profile, lightweight base foam, blue XXSoft FRG foam topper., CoolCore® cushion cover	 High-density foam base, more contour to sacral/coccyx relief well, Stealth Tek, 2-way stretch breathable cover (fluid resistant top and side panels)
WEIGHT CAPACITY	14-15" W: 175 lbs. 16-21" W: 350 lbs. 22+" W: 450 lbs.	10-15" W: 175 lbs. 16-21" W: 350 lbs. 22" W: 450 lbs.	10-15" W: 175 lbs. 16-20" W: 350 lbs. 22+" W: 500 lbs.	14" W: 175 lbs. Foam & Gel 16-20" W: 350 lbs. Foam & Gel 22+" W: 650 lbs. Foam; 500 lbs. Gel	14-15" W: 175 lbs. 16-21" W: 350 lbs. 22" W: 500 lbs.	14-15" W: 175 lbs. 16-21" W: 350 lbs. 22+" W: 400 lbs.	14-15" W: 175 lbs. 16-21" W: 350 lbs. 22+" W: 400 lbs.	16-21" W: 350 lbs. 22+" W: 400 lbs.
WEIGHT⁴	3.9 lbs.	5.5 lbs.	4.9 lbs.	Foam: 3.75 lbs. Gel: 4.3 lbs.	3.55 lbs.	4.3 lbs.	2.5 lbs.	2.2 lbs.
HEIGHT³ (AT HIGHEST POINT OF CONTOUR)	4.9"	3.7"	3"	3"	3.7"	2.9"	2.9"	2.9"
DIMENSIONS^{2,3}	14x14" up to 24x24" Odd sizes: Widths 15, 17, 19, 21" Depths 15, 17, 19, 21" SOF available	10x10" up to 22x24" Odd sizes: Widths 15, 17, 19, 21" Depths 15, 17, 19, 21" SOF available	10x10" up to 24x24" Odd sizes: Widths 15, 17, 19, 21" Depths 15, 17, 19" SOF available	Foam: 14x18" up to 32x28" Gel: 14x18" up to 24x26" Odd sizes: Widths 17, 19" Depths 17, 19, 21" SOF available	14x13" up to 22x22" Odd sizes: Widths 15, 17, 19" Depths 15, 17, 19" SOF available	14x12" up to 24x24" Odd sizes: Widths 15, 17, 19, 21" Depths 15, 17, 19, 21, 23" No SOF	14x12" up to 24x24" Odd sizes: Widths 15, 17, 19, 21" Depths 15, 17, 19, 21, 23" SOF available	16x16" up to 24x24" Odd sizes in 15-21" W No SOF
MSRP	\$515 - \$575	Essence Basic: \$485 - \$435 Essence Plus: \$725 - \$775	\$580 - \$645	Spectrum® Foam: \$580 - \$630 Spectrum® Gel: \$580 - \$650	\$585 - \$730	\$265 - \$530	\$455 - \$765	\$310 - \$370
ADP	(sizes up to 20"x20") SE0001524 (sizes > 20"x20") SE0001525	Essence Basic: (sizes up to 20" X 20") SE0001481 (sizes > 20" X 20") SE0001482 Essence Plus: (sizes up to 20" X 20") SE0001487 (sizes > 20" X 20") SE0001488	(sizes up to 20"x20") SE0001526 (sizes > 20"x20") SE0001527	Spectrum® Foam: (sizes up to 20" X 20") SE0001528 (sizes > 20" X 20") SE0001529 Spectrum® Gel: (sizes up to 20" X 20") SE0001530 (sizes > 20" X 20") SE0001531	(sizes 16"-20" wide) SE0001622 (sizes > 20" wide) SE0001623	(sizes up to 15" wide) SE0001619 (sizes 16-19" wide) SE0001620 (sizes 20-24" wide) SE0001621	(sizes 10-15" wide) SE0001651 (sizes 16-20" wide) SE0001652 (sizes 21-24" wide) SE0001653	(sizes up to 20"x20") SE0001522 (sizes > 20"x20") SE0001523

1) CoolCore® is standard on all cushions except the Simplicity G.
 2) For a complete list of available sizes refer to the cushion order form. For special order sizes contact sales for assistance. Special order sizes available for all cushions except Glacial SP and the Simplicity G.
 3) Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.
 4) Each weight corresponds to an 18"x18" cushion.