



STEALTH POSITIONING

Stability Is a Snap



stealthproducts.com



Stealth's Positioning line of pelvic belts, vests, chest straps and foot positioners are intended to increase trunk stability, aid in properly aligning the pelvis, and provide correct posture in a wheelchair. With Stealth Products, stability is a snap!



Breathable Material

Constructed from a unique open-cell material that allows air to circulate and moisture to escape. This contributes to reduced heat buildup and a cooler, drier, and more comfortable user experience.



Dynamic Support

Dynamic Support provides postural control, a degree of trunk movement, and a softer fit on the body. The user has movement within a limited range encouraging return to their neutral position.



Offset Belt Anchors

A slight offset allows for easy adjustment without loosening the anchor bolts.



Structured Support

The structured material is more resistant to high tone which helps maintain correct posture.



Swivel Buckles

Swivel buckles on the lower attaching straps allow self-adjustment, supporting the correct angle of attachment.



Full Length Attaching Straps

Provide for proper positioning of the top side release buckles, ensuring that the posture support maintains adjustment.




Posture Support

- Breathable Material
- Dynamic & Structured Support
- Classic & Contoured Styles
- Swivel Buckles
- Full Range of Sizes
- Offset Belt Anchor
- Large Neck Cavity
- Full Length Attaching Straps

Classic Style

The Classic Style offers a large padded surface for the thoracic area.




	A	B	C	D	E	F	G	H
XS	7.5	12.5	2	3.5	7.5	3.5	3	2
S	8.5	15	2	4.5	9	4	4	2
M	9.5	19	2	5.5	11	5	5	2.75
L	10.5	22.5	2	6.5	13	6	6	3.5
XL	11.5	25.5	2.5	6.5	14	7	7	4.5

	XS	S	M	L	XL
Dynamic	IFPS1101	IFPS1102	IFPS1103	IFPS1104	IFPS1105
Structured	IFPS1201	IFPS1202	IFPS1203	IFPS1204	IFPS1205

Contoured Style

The Contoured Style has a narrowed section over the sternum that allows a better fit for females, G-tube users or those requiring a smaller support surface.



	A	B	C	D	E	F	G	H
XS	7.5	12.5	2	3.5	7.5	1	3.5	1.25
S	8.5	15	2	4.5	9	1	4.25	1.75
M	9.5	18.75	2	5.5	11	1	5.25	2.5
L	10.5	22.5	2	6.5	13	1.25	6	3.5
XL	11.5	25.5	2.5	6.5	14	1.5	7	4.5

	XS	S	M	L	XL
Dynamic	IFPS2101	IFPS2102	IFPS2103	IFPS2104	IFPS2105
Structured	IFPS2201	IFPS2202	IFPS2203	IFPS2204	IFPS2205



CLASSIC STYLE VINYL FACE (DYNAMIC)

IFPS1401	Pull Down, Extra Small
IFPS1402	Pull Down, Small
IFPS1403	Pull Down, Medium
IFPS1404	Pull Down, Large
IFPS1405	Pull Down, Extra Large

CLASSIC STYLE VINYL FACE (STRUCTURED)

IFPS1301	Pull Down, Extra Small
IFPS1302	Pull Down, Small
IFPS1303	Pull Down, Medium
IFPS1304	Pull Down, Large
IFPS1305	Pull Down, Extra Large



CONTOURED STYLE VINYL FACE (DYNAMIC)

IFPS2401	Pull Down, Extra Small
IFPS2402	Pull Down, Small
IFPS2403	Pull Down, Medium
IFPS2404	Pull Down, Large
IFPS2405	Pull Down, Extra Large

CONTOURED STYLE VINYL FACE (STRUCTURED)

IFPS2301	Pull Down, Extra Small
IFPS2302	Pull Down, Small
IFPS2303	Pull Down, Medium
IFPS2304	Pull Down, Large
IFPS2305	Pull Down, Extra Large



CONTOURED STYLE ZIPPER FACE (DYNAMIC)

IFPS3101	Pull Down, Extra Small
IFPS3102	Pull Down, Small
IFPS3103	Pull Down, Medium
IFPS3104	Pull Down, Large
IFPS3105	Pull Down, Extra Large

CLASSIC STYLE ZIPPER FACE (STRUCTURED)

IFPS3201	Pull Down, Extra Small
IFPS3202	Pull Down, Small
IFPS3203	Pull Down, Medium
IFPS3204	Pull Down, Large
IFPS3205	Pull Down, Extra Large

All Sizes Are Inches



BACKPACK STYLE STRAP, PAIR (STRUCTURED)

IFPB4201 Pull Down, Extra Small

IFPB4202 Pull Down, Small

IFPB4203 Pull Down, Medium

IFPB4204 Pull Down, Large

H-STYLE HARNESS (STRUCTURED)



IFHH5201 Pull Down, Extra Small

IFHH5202 Pull Down, Small

IFHH5203 Pull Down, Medium

IFHH5204 Pull Down, Large

Chest Straps (All chest straps are dynamic)




		A	B	C
Belt Anchor	XS	14	2.8	3.5
	S	16	2.8	3.75
	M	18	3.1	4
	L	21	3.1	4.5
	XL	23	3.25	4.75

		A	B
Velcro	XS	7	2
	S	8.5	3
	M	10.25	4
	L	11.5	5
	XL	13	6

	XS	S	M	L	XL
Belt Anchor	IFCS4111	IFCS4112	IFCS4113	IFCS4114	IFCS4115
Velcro	IFCS4121	IFCS4122	IFCS4123	IFCS4124	IFCS4125

Y Straps



	A	B	C
S	19	8	2
M	22	9	2.25
L	25	11	2.5


	S	M	L
Dynamic	IFYS3102	IFYS3103	IFYS3104
Structured	IFYS3202	IFYS3203	IFYS3204



Pelvic Belts & Pads


- Help maintain a neutral pelvic position, which directly affects seated body alignment and distributes pressure appropriately.
- Increase postural stability, decreasing the effort needed to sit upright and facilitating upper extremity movement.
- Contribute to proper head positioning and neck alignment.
- Reduce negative postural tendencies that can lead to orthopedic asymmetries and other secondary medical complications.

2 Point Pad Sizes




	Part Number	A	B	C
XS	IFPB101	4.25	1.75	1.25
S	IFPB102	5.5	1.75	1.25
M	IFPB103	7.5	2.5	1.75
L	IFPB104	9.5	2.75	2.25
XL	IFPB105	11.5	2.75	2.25

4 Point Pad Sizes




	Part Number	A	B	C
XS		4.25	1.75	3.25
S		5.5	2.25	4
M		7.5	2.5	4
L		9.5	2.25	4

Belt Sizes



	A	B	C	Buckle	Pad Sizes
Single Pull	32	27	1	PB	NP / XS
Single Pull	32	27	1	SR	NP
Dual Pull	24	24	1	PB	XS / S
Single Pull	32	27	1.5	PB / SR	NP / M / L
Dual Pull	24	24	1.5	PB / SR	M / L
Single Pull	32	27	2	PB	NP / L / XL
Single Pull	32	27	2	SR	L / XL
Dual Pull	24	24	2	PB / SR	L / XL



	A	B	C	D	Buckle	Sizes
	32	27	20	1	PB	XS
	32	27	20	1.5	PB	S / M / L

Part Numbers



2 POINT NON-PADDED BELTS (PUSH BUTTON)

- IFPB5110 Push Button, 1 in, Single Pull, Non-Padded
- IFPB5210 Push Button, 1.5 in, Single Pull, Non-Padded
- IFPB5310 Push Button, 2 in, Single Pull, Non-Padded
- IFPB5350 Pelvic Belt, Push Button, 2 in, Retractable, Non-Padded

2 POINT NON-PADDED (SIDE RELEASE)

- IFPB6110 Side Release, 1 in, Single Pull, Non-Padded (shown above)
- IFPB6210 Side Release, 1.5 in, Single Pull, Non-Padded



2 POINT PADDED BELTS (DUAL PULL)

- IFPB5121 Push Button, 1 in, Dual Pull, Extra Small Pads (4.5 in)
- IFPB5122 Push Button, 1 in, Dual Pull, Small Pads (5.5 in)

1 1/2 inch belts

- IFPB5223 Push Button, 1.5 in, Dual Pull, Medium Pads (7.5 in)
- IFPB6223 Side Release, 1.5 in, Dual Pull, Medium Pads (7.5 in)
- IFPB5224 Push Button, 1.5 in, Dual Pull, Large Pads (9.5 in)
- IFPB6224 Side Release, 1.5 in, Dual Pull, Large Pads (9.5 in)

2 inch belts

- IFPB5324 Push Button, 2 in, Dual Pull, Large Pads (9.5 in)
- IFPB6324 Side Release, 2 in, Dual Pull, Large Pads (9.5 in)
- IFPB5325 Push Button, 2 in, Dual Pull, Extra Large Pads (11.5 in)
- IFPB6325 Side Release, 2 in, Dual Pull, Extra Large Pads (11.5 in)



2 POINT PADDED BELTS (SINGLE PULL)

- IFPB5111 Push Button, 1 in, Single Pull, Extra Small Pads (4.5 in)
- IFPB5112 Push Button, 1 in, Single Pull, Small Pads (5.5 in)

1 1/2 inch belts

- IFPB5213 Push Button, 1.5 in, Single Pull, Medium Pads (7.5 in) (shown above)
- IFPB6213 Side Release, 1.5 in, Single Pull, Medium Pads (7.5 in)
- IFPB5214 Push Button, 1.5 in, Single Pull, Large Pads (9.5 in)
- IFPB6214 Side Release, 1.5 in, Single Pull, Large Pads (9.5 in)

2 inch belts

- IFPB5314 Push Button, 2 in, Single Pull, Large Pads (9.5 in)
- IFPB6314 Side Release, 2 in, Single Pull, Large Pads (9.5 in)
- IFPB5315 Push Button, 2 in, Single Pull, Extra Large Pads (11.5 in)
- IFPB6315 Side Release, 2 in, Single Pull, Extra Large Pads (11.5 in)

2 POINT PADDED BELTS (REAR PULL)

- IFPB5131 Push Button, 1 in, Rear Pull, Extra Small Pads (4.5 in) (E0978)
- IFPB5132 Push Button, 1 in, Rear Pull, Small Pads (5.5 in) (E0978)

1 1/2 inch belts

- IFPB5233 Push Button, 1.5 in, Rear Pull, Medium Pads (7.5 in) (E0978)
- IFPB5234 Push Button, 1.5 in, Rear Pull, Large Pads (9.5 in) (E0978)
- IFPB6233 Side Release, 1.5 in, Rear Pull, Medium Pads (7.5 in) (E0978)
- IFPB6234 Side Release, 1.5 in, Rear Pull, Large Pads (9.5 in) (E0978)

2 inch belts

- IFPB5334 Push Button, 2 in, Rear Pull, Large Pads (9.5 in) (E0978)
- IFPB5335 Push Button, 2 in, Rear Pull, Extra Large Pads (11.5 in) (E0978)
- IFPB6334 Side Release, 2 in, Rear Pull, Large Pads (9.5 in) (E0978)
- IFPB6335 Side Release, 2 in, Rear Pull, Extra Large Pads (11.5 in) (E0978)



4 POINT BELTS (V-CONTOUR STYLE)

- IFPB5141 Push Button, 1 in, 4-Point, Extra Small Pads (4.5 in)

1 1/2 inch belts

- IFPB5242 Push Button, 1.5 in, 4-Point, Small Pads (5.5 in)
- IFPB5243 Push Button, 1.5 in, 4-Point, Medium Pads (7.5 in)
- IFPB5244 Push Button, 1.5 in, 4-Point, Large Pads (9.5 in)



4 POINT BELTS (LINEAR STYLE)

- IFPB5145 Push Button, 1 in, 4-Point, Extra Small Pads (4.5 in)

1 1/2 inch belts

- IFPB5246 Push Button, 1.5 in, 4-Point, Small Pads (5.5 in)
- IFPB5247 Push Button, 1.5 in, 4-Point, Medium Pads (7.5 in)
- IFPB5248 Push Button, 1.5 in, 4-Point, Large Pads (9.5 in)




RETRACTABLE SEATBELT

- IFRB5350 Modification, Pelvic Belt, Push Button, Retractable 2" Web, Non-Padded

Foot Positioning

Ankle Support Sizes



	A (Circumference)	B (Width)	C (Length)
XS	6-8		
S	7-8.5	4.75"	7.75"
M	8-11	5.25"	9.25"
L	9-12	5.75"	10.25"
XL	12.5-15		



ANKLE POSITIONING SYSTEMS

IFFP1110	Dynamic 3/4 in Webbing, Extra Small
IFFP1111	Dynamic 3/4 in Webbing, Small
IFFP1122	Dynamic 1 in Webbing, Medium
IFFP1123	Dynamic in Webbing, Large
IFFP1124	Dynamic 1 in Webbing, Extra Large



FULL FOOT POSITIONING

IFFP2110	Dynamic 3/4 in Webbing, Extra Small
IFFP2111	Dynamic 3/4 in Webbing, Small
IFFP2122	Dynamic 1 in Webbing, Medium
IFFP2123	Dynamic 1 in Webbing, Large
IFFP2124	Dynamic 1 in Webbing, Extra Large




FIXED ANKLE SUPPORTS

IFFP3210	Structured 3/4 in Webbing, Extra Small
IFFP3211	Structured 3/4 in Webbing, Small
IFFP3212	Structured 1 in Webbing, Medium
IFFP3213	Structured 1 in Webbing, Large
IFFP3214	Structured 1 in Webbing, Extra Large

FIXED ANKLE SUPPORTS (VELCRO)

IFFP3310	Fixed Ankle Support, Structured Velcro, 3/4 in Webbing, Extra Small
IFFP3311	Fixed Ankle Support, Structured Velcro, 3/4 in Webbing, Small
IFFP3312	Fixed Ankle Support, Structured Velcro, 1 in Webbing, Medium
IFFP3313	Fixed Ankle Support, Structured Velcro, 1 in Webbing, Large
IFFP3314	Fixed Ankle Support, Structured Velcro, 1 in Webbing, Extra Large

Shoe Holder Sizes



	A	B
S	5	3
M	6.25	3.5
L	7.75	4.25



SHOE HOLDERS

IFSH4211	Shoe Holder, Padded Straps 3/4 in Webbing, S
IFSH4212	Shoe Holder, Padded Straps 3/4 in Webbing, M
IFSH4213	Shoe Holder, Padded Straps 3/4 in Webbing, L
IFTS4211	Padded Toe Strap Only, S
IFTS4212	Padded Toe Strap Only, M
IFTS4213	Padded Toe Strap Only, L
IFSH4001	Shoe Holder Pad, Velcro, S
IFSH4002	Shoe Holder Pad, Velcro, M
IFSH4003	Shoe Holder Pad, Velcro, L

Leg Support

LEG HARNESS

IFLH1111	Leg Harness, Dynamic, Small, 1.5 in Webbing
IFLH1112	Leg Harness, Dynamic, Medium, 1.5 in Webbing
IFLH1123	Leg Harness, Dynamic, Large, 2 in Webbing

LEG HARNESS



CALF STRAP



CALF PANEL



Calf Support

CALF STRAPS (STRUCTURED)

IFCF1511	Calf Strap, Structured, Small
IFCF1512	Calf Strap, Structured, Medium
IFCF1513	Calf Strap, Structured, Large

CALF STRAPS (DYNAMIC)

IFCF1611	Calf Strap, Dynamic, Small
IFCF1612	Calf Strap, Dynamic, Medium
IFCF1613	Calf Strap, Dynamic, Large

CALF PANELS (STRUCTURED)

IFCF2511	Calf Panel, Structured, Small
IFCF2512	Calf Panel, Structured, Medium
IFCF2513	Calf Panel, Structured, Large
IFCF2514	Calf Panel, Structured, Extra Large

CALF PANELS (DYNAMIC)

IFCF2611	Calf Panel, Dynamic, Small
IFCF2612	Calf Panel, Dynamic, Medium
IFCF2613	Calf Panel, Dynamic, Large
IFCF2614	Calf Panel, Dynamic, Extra Large

All Sizes Are Inches

Components

**Padded
Toe Strap**



IFFP0110	Dynamic Ankle Support REPLACEMENT , XS
IFFP0111	Dynamic Ankle Support REPLACEMENT , S
IFFP0122	Dynamic Ankle Support REPLACEMENT , M
IFFP0123	Dynamic Ankle Support REPLACEMENT , L
IFFP0124	Dynamic Ankle Support REPLACEMENT , XL
IFTS2111	Padded Toe Strap, Foot Positioning System, S
IFTS2122	Padded Toe Strap, Foot Positioning System, M
IFTS2123	Padded Toe Strap, Foot Positioning System, L
IFFP100-01	Foot & Ankle Base Plate REPLACEMENT , S
IFFP100-02	Foot & Ankle Base Plate REPLACEMENT , M
IFFP100-03	Foot & Ankle Base Plate REPLACEMENT , L
IFTS0322	Toe Strap, Hook and Loop 1" x 12", One Size

Accessories



**Buckle
Cover**



**Cane
Clamp**

IFMDBC10	Modification, Buckle Cover Kit. 1 in Buckle
IFMDBC15	Modification, Buckle Cover Kit. 1.5 in Buckle
IFMDCM	Modification, Cam Buckle Mount Kit
IFMDCC	Modification, Cane Clamp Kit



Distributed by: Pride Mobility Products Co.
5096 South Service Rd. Beamsville, Ontario L0R 1B3
<https://www.pridemobility.ca/>
[stealthproducts.com](https://www.pridemobility.ca/)