

# TRU BALANCE<sup>®</sup>

## FLEX



Introducing **TRU-Balance<sup>®</sup> Flex Seating**,  
a new low-cost static seat option  
available on the Edge 3 Stretto<sup>™</sup>, Edge 3,  
Q6 Edge<sup>®</sup> 2.0, Q4 and the Rival<sup>®</sup>.

TRU-Balance<sup>®</sup> Flex Seating  
shown on Q6 Edge<sup>®</sup> 2.0 Power Base  
with Q-Logic 2 Controller



**QUANTUM<sup>®</sup>**  
#1 FOR REHAB POWER

Beamsville, Ontario

| 888-570-1113

| [quantumrehab.ca](http://quantumrehab.ca)



## SPECIFICATIONS

<b>Weight Capacity:</b>	300 lbs.
<b>Seat Dimensions:</b>	<b>Width:</b> 15-22" <b>Depth:</b> 15-20"
<b>Manual Recline:</b>	-5° - 30°
<b>Flip-Forward Back Cane:</b>	Angle can be set from 85° to 115° (in 5° increments)
<b>Back Options:</b>	SitRite Premier
<b>Seat-to-Floor Height:</b>	<b>Q6 Edge® 2/Q4:</b> 16.5" <b>Rival:</b> (3 Settings) 16.5", 17.3", 18.1"
<b>Arms:</b>	Flex drop-in "T" arm (Desk or Full Length)
<b>Arm Height Ranges:</b>	9"-14" adjustable in 1" increments

## FEATURES

- Easily adjustable, this versatile design can grow with your client
- Same ease of width and depth adjustment as TRU-Balance® 3
- No-charge, TRU-Balance® Flex, height-adjustable, single-post armrests with clothing guard standard
- Accepts all rail-mounted TRU-Balance 3 accessories
- Angle adjustable fold down back canes
- Single seat pan covers all sizes
- Height and angle adjustable seat towers

## OPTIONAL COMPONENTS

- TRU-Balance® 3 style back canes  
Options include:
  - standard angle adjustable hatchet assembly
  - quick fold-down cane assembly
  - Endo and Meso offset brackets
  - optional height from 14"-25"
- TRU-Balance® Flex manual recline system
- TRU-Balance Flex back assembly with SitRite Premier back cushion
- Accepts aftermarket back rests
- TRU-Balance single-post, dual-post and pediatric arm rests
- Cane mounted flip-back arms – available on TRU-Balance canes
- TRU-Balance #8 swing away legrests
- Elevating legrests
- Sling upholstery – available on TRU-Balance canes
- TRU-Comfort® 2 back system
- Synergy® Backs – available on TRU-Balance canes

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.